

# Optimizing Preventive Cardiology: Harnessing AI For Early Detection Of Coronary Artery Disease

Yotam Kimmel, B.Sc.<sup>1</sup>, Alexis Kurek, P.hD.<sup>2</sup>, Orit Wimpfheimer, MD<sup>1</sup>, Oren Shalem, B.Sc.<sup>1</sup>, David Langholz, MD<sup>2</sup>

<sup>1</sup> **Nanox.AI**, Petach Tikva, Israel | <sup>2</sup> **Corewell Health**, Grand Rapids, MI, USA.

## Introduction

- Coronary artery disease (CAD)** accounts for 50% of all cardiovascular related deaths and remains a leading cause of death worldwide. Whilst cardiac gated CTs minimize image noise caused by cardiac motion, there are an estimated 20 times as many non-cardiac-gated CTs of the chest performed in the USA each year.
- The HealthCCSng (NanoxAI) device is an **artificial intelligence software** to evaluate calcified plaques in the coronary arteries on non-ECG gated chest CTs, which may present risk for coronary artery disease. The detection capabilities of AI make it a valuable tool in population health management, as organizations work to shift their services to early identification and intervention of chronic disease.
- The cascade of data from identification to care delivery represents a challenge to many health care systems. Corewell Health **has successfully deployed a system integrating the AI data** into the EMR for the early detection and intervention of CAD.

## Methods

- The device was deployed in December 2022 at Corewell Health for all CT scans meeting inclusion criteria for 14months. Studies with high or medium coronary artery calcium (CAC) (CAC>99 AU) were sent for radiologist review.
- Text related to the finding was embedded discretely into the radiology report. This text was then filed into the EMR, triggering a chart query to determine if the patient had an existing diagnosis of CAD.
- The cascade of data from identification to care delivery represents a challenge to many health care systems. Corewell Health **has successfully deployed a system integrating the AI data** into the EMR for the early detection and intervention of CAD.



## Results

- Data was collected between 2/1/23 and 1/31/24. **32,650 scans** were read by the algorithm, with **8,481 returning medium or high CAC findings**; **67% were reported** by the radiologist. Of these, **64% were previously undiagnosed** with CAD and **79% were not on statin therapy**. Those undiagnosed had an average age of 70.8, an average End of Life score of 11.2, and an average ASCVD risk of 19.1.
- History of tobacco use and oncological concerns were the most frequent ordering diagnoses. Following identification, **14.3% of patients were formally diagnosed with CAD**, **27% were started on statins**, and **12% were referred to cardiology**.


## Conclusions

- This observational study demonstrated that the AI algorithm is useful in **identifying patients with undiagnosed CAD** and that AI integrated into the EMR **positively impacts population health goals**.

### Sample output from the NanoxAI HealthCCSng



### Diagnosed patients



### Patient journey CAD diagnosis

- 812 Patients** Formally diagnosed with CAD following NanoxAI results
- 1534 Patients** Started on statins following NanoxAI results
- 682 Patients** Were referred to cardiology following NanoxAI results