

# Leveraging AI for Opportunistic Screening: Identifying Coronary Artery Calcification on Non-ECG Gated Lung Cancer Screening Chest CT

Kevin Pham DO, Akash Dadlani MD, Ava Wexler MD, Devon Klein MD MPH



Atlantic  
Health System

# Disclosures

Nanox AI provided software free of charge for the purposes of conducting this study.

# Background

- Atherosclerotic cardiovascular disease (ASCVD) is prevalent in the general population and remains one of the leading causes of morbidity and mortality in the United States
- Coronary artery calcium (CAC) is a highly specific feature of coronary atherosclerosis which has been shown to predict mortality
- With advances in computed tomography, CAC scoring has become a widely available way to assess ASCVD risk for major cardiovascular events, particularly for those who are asymptomatic
- CAC scoring has shown to be cost effective as its detection allows for treatment through primary prevention strategies

# Purpose

- Asymptomatic adults with intermediate to high ASCVD risk are typically screened for CAC with dedicated ECG-gated chest CT examinations
- CAC can be detected incidentally on non-ECG gated chest CT studies performed for other reasons including trauma, dyspnea, and lung cancer screening
- Evaluation of the presence of CAC and its severity on all non-contrast chest CTs is recommended by the Society of Cardiovascular Computed Tomography and Society of Thoracic Radiology based on their 2016 guidelines
- Using commercially available artificial intelligence (AI) solutions, we can identify patients with moderate or severe coronary artery calcification with high specificity
- Early detection of CAC can be clinically impactful for patients without known coronary artery disease (CAD) who would otherwise not be screened

# Current State of Risk Assessment on Non-Dedicated Exams



Presence of coronary calcification may not be reported



When reported, there is no standard assessment

Binary reporting of presence or absence

Subjective “eyeball” estimation of mild-moderate-severe

## Examinations Across Atlantic Health System



Approximately 10,000 outpatient non-contrast chest CTs were performed within our health system in 2023



Of these, 3200 were low-dose chest CTs performed for lung cancer screening



Utilizing Nanox AI, we can identify and report on the presence and severity of CAC more consistently on these non-targeted studies

# Nanox AI Cardiac Solution


Analyzes non-contrast, non-ECG gated chest CTs

Automatically measures calcified plaque

- Quantifies by Agatston score
- Agatston risk stratification based on total calcium score
  - Zero
  - Minimal (1-10)
  - Mild (11-100)
  - Moderate (101-400)
  - Severe (>400)

Outputs key images for verification

# Nanox AI Workflow





# Nanox AI Examination Criteria

## Inclusion Criteria

- Modality: CT
- Gender: All
- Orientation: Axial
- Age: 30 years old and older
- Slice thickness 0.1 mm – 3.1 mm (inclusive)
- Slice interval 0.1 mm – 3.1 mm (inclusive)
- Non-contrast studies only
- Any CT that includes the heart

## Exclusion Criteria

- Any CT not including the entire heart
- ECG-gated CTs
- Number of slices: <20
- Series: CT Attenuation Correction, PET

# Nanox AI Sample Image

**Study Information**



Patient ID: ANONAL80GI111    Patient Name: , Anonymous M...    Gender: M    Study Date: 07/09/2023  
DOB: 29/01/1950    Age: 73    Accession #: RAL80GI1T0

**Feedback**

Agree     Disagree    + Add comment

**Result:** Medium

**Detected Calcium**






Not for diagnostic use    2/5

Index bar indicates CAC quantification

Segmentation of coronary artery indicated in yellow

# Objective Coronary Calcium Scoring



# Two Phase Pilot at Atlantic Health System

- Phase 1
  - Validation
    - Does the AI work with our population and technology?
    - How user friendly is it?
    - Do our results match published data?
  
- Phase 2
  - Experimental Phase
    - Is it reproducible if we apply it to a subset of our population?
    - What can we learn?
      - ♦ What is the incidence of undocumented disease?
      - ♦ Can we derive any demographic relationships?
      - ♦ How do we leverage our findings to improve the health of our patients?

## Phase 1: Validation

- 497 adult non-contrast, non-ECG gated chest CTs were retrospectively submitted to Nanox AI for analysis
- Of these, 382 (77%) studies were successfully analyzed by the AI algorithm
- In 41.6% of cases (159/382), moderate to severe disease was discovered
  - 30.1% with CAC Score > 400 (Severe)
  - 11.5% with CAC Score of 100-400 (Moderate)

## Phase 2: Experimental

- 492 lung cancer screening studies were submitted for analysis by Nanox AI
- 87.6% (431/492) were successfully analyzed
- 39.7% (171/431) of patients were found to have moderate to severe disease
  - 23.4% (101) with CAC Score > 400 (Severe)
  - 16.2% (70) with CAC Score of 100-400 (Moderate)

# Study Population Demographics

- Mean patient age = 67 (SD = 5.8)
- Age range = 52 – 78

Gender	Frequency	Percentage
Male	120	70
Female	51	30
Total	171	100

# Chart Review

- AI detected moderate to severe CAC burden in 39.7% (171) patients undergoing lung cancer screening with low-dose chest CT
- Upon chart review, 31.0% (53) of these patients did not have a known diagnosis of coronary artery disease at the time of imaging
- Known/existing coronary artery disease was defined as a documented diagnosis of CAD, documented statin therapy, and/or a history of coronary artery stenting

At the time of imaging, did the patient have known/existing CAD?	Frequency	Percentage
Yes	118	69
No	53	31
Total	171	100



## Lung-RADS Score vs Known/Existing CAD at Imaging

Lung RADS	Known/Existing CAD at Imaging?		Total
	No	Yes	
3	3 30%	7 70%	10 100%
3S	2 67%	1 33%	3 100%
4A	1 17%	5 83%	6 100%
4B	2 100%	0 0%	2 100%
4X	1 100%	0 0%	1 100%
Total	9 41%	13 59%	22 100%

# Correlation Between Lung-RADS Score and Moderate-Severe CAC

Lung RADS	Frequency	Percentage
1	31	18.1
1S	4	2.3
2	97	56.7
2S	17	9.9
3	10	5.8
3S	3	1.8
4A	6	3.5
4B	2	1.2
4X	1	0.6
Total	171	100

A Lung-RADS score of 3 or higher is considered significant

13% (22/171) of patients with moderate to severe CAC also had significant Lung-RADS findings

# Correlation Between Lung-RADS Score and Moderate-Severe CAC

Lung RADS	AI Result of Disease		Total
	Moderate	Severe	
Less than 3	61 41%	88 59%	149 100%
Greater than or equal to 3	9 41%	13 59%	22 100%
Total	70 41%	101 59%	171 100%

No significant relationship was seen between Lung-RADS Score and the presence of moderate to severe CAC

# Relationship Between Gender and AI Results

Gender	AI Result of Disease		Total
	Moderate	Severe	
Female	26 51%	25 49%	51 100%
Male	44 37%	76 63%	120 100%
Total	70 41%	101 59%	171 100%

Male patients undergoing lung cancer screening were found to have a higher incidence of severe CAC burden compared to female patients

# Conclusion

- Using a commercially available AI solution, we were able to successfully identify the presence of CAD in a significant number of patients within our health system who did not previously have a known or existing CAD
- Moderate to severe CAC burden was detected by AI in 41.6% and 39.7% of cases following each phase of the study
- 31.0% of patients undergoing lung cancer screening with low-dose chest CT were found to have moderate to severe CAC and did not have a known or existing diagnosis of CAD at the time of imaging
- No significant relationship was demonstrated between moderate to severe CAC burden and significant Lung-RADS findings
- Utilizing AI to detect incidental CAC in patients undergoing non-dedicated CT exams can lead to early detection of clinically significant CAD
  - This can lead to early intervention with primary prevention strategies which would ultimately improve long-term outcomes and decrease healthcare costs

# References

- Agatston AS, Janowitz WR, Hildner FJ, Zusmer NR, Viamonte M, Detrano R. Quantification of coronary artery calcium using ultrafast computed tomography. *Journal of the American College of Cardiology*. 1990;15(4):827-832. doi:10.1016/0735-1097(90)90282-t
- Greenland P, Blaha MJ, Budoff MJ, Erbel R, Watson KE. Coronary Calcium Score and Cardiovascular Risk. *J Am Coll Cardiol*. 2018 Jul 24;72(4):434-447. doi: 10.1016/j.jacc.2018.05.027
- Hecht HS, Cronin P, Blaha MJ, et al. 2016 SCCT/STR guidelines for coronary artery calcium scoring of noncontrast noncardiac chest CT scans: A report of the Society of Cardiovascular Computed Tomography and Society of Thoracic Radiology. *Journal of Cardiovascular Computed Tomography*. 2017;11(1):74-84. doi:10.1016/j.jcct.2016.11.003
- Kerndt CC, Chopra R, Weber P, Rechenberg A, Summers D, Boyden T, Langholz D. Using Artificial Intelligence to Semi-Quantitate Coronary Calcium as an 'Incidentaloma' on Non-Gated, Non-Contrast CT Scans, A Single-Center Descriptive Study in West Michigan. *Spartan Med Res J*. 2023 Dec 5;8(1):89132. doi: 10.51894/001c.89132
- Obisesan OH, Osei AD, Uddin SMI, Dzaye O, Blaha MJ. An Update on Coronary Artery Calcium Interpretation at Chest and Cardiac CT. *Radiol Cardiothorac Imaging*. 2021 Feb 25;3(1):e200484. doi: 10.1148/ryct.2021200484
- Park C, Lee BC, Jeong WG, Park W-J, Jin GY, Kim Y-H. Coronary artery calcification on low-dose lung cancer screening CT in South Korea: Visual and Artificial Intelligence–Based Assessment and association with Cardiovascular events. *American Journal of Roentgenology*. 2024;222(5). doi:10.2214/ajr.24.30852
- Sandhu AT, Rodriguez F, Ngo S, et al. Incidental coronary artery calcium: Opportunistic screening of previous nongated chest computed tomography scans to improve statin rates (notify-1 project). *Circulation*. 2023;147(9):703-714. doi:10.1161/circulationaha.122.062746
- Shemesh J, Henschke CI, Shaham D, Yip R, Farooqi AO, Cham MD, McCauley DI, Chen M, Smith JP, Libby DM, Pasmantier MW, Yankelevitz DF. Ordinal scoring of coronary artery calcifications on low-dose CT scans of the chest is predictive of death from cardiovascular disease. *Radiology*. 2010 Nov;257(2):541-8. doi: 10.1148/radiol.10100383
- Xie X, Zhao Y, de Bock GH, de Jong PA, Mali WP, Oudkerk M, Vliegenthart R. Validation and prognosis of coronary artery calcium scoring in nontriggered thoracic computed tomography: systematic review and meta-analysis. *Circ Cardiovasc Imaging*. 2013 Jul;6(4):514-21. doi: 10.1161/CIRCIMAGING.113.000092